



REFERENCE GUIDE

Side 1	<p>Reference Guide Repatha® (ri-PATH-a) (evolocumab) injection, for subcutaneous use Single-Dose Prefilled SureClick® Autoinjector 140 mg/mL</p>
---------------	--

Read all instructions in carton before use.


Questions? Call 1-844-REPATHA (1-844-737-2842)

Guide to parts	
Gray start button Expiration date Window Medicine Orange cap on	
<p>Wait at least 30 minutes for the autoinjector to come to room temperature before giving your injection.</p>	
<p>Prepare and clean your injection site.</p>	
1	

Pull the orange cap straight off only when you are ready to inject.


Do not leave the orange cap off for **more than 5 minutes**. This can dry out the medicine.

2










Orange cap

5 minutes



Turn over to continue

<p>Side 2</p>	<p>Reference Guide Repatha® (ri-PATH-a) (evolocumab) injection, for subcutaneous use Single-Dose Prefilled SureClick® Autoinjector 140 mg/mL</p>	<p>AMGEN[®] Manufactured by: Amgen Inc. One Amgen Center Drive Thousand Oaks, CA 91320-1799 © 2019-2022 Amgen Inc. All rights reserved.  This printed material is recyclable. <part number> Revised: 8/2022 v4</p>	<p>Read other side first</p>	
<p>Stretch or pinch your injection site to create a firm surface.</p>		<p>Put the yellow safety guard (needle inside) on your skin at 90 degrees.</p>	<p>Firmly push down the autoinjector onto your skin until it stops moving.</p>	<p>Keep pushing down on your skin for about 15 seconds.</p>
<p>3</p> <p>Stretch method:</p>  <p>or</p> <p>Pinch method:</p> 	<p>4</p> <p>Yellow safety guard (needle inside)</p> 	<p>5</p> <p>When you are ready to inject, press the gray start button. You will hear a click.</p> <p>"click"</p> 	<p>6</p> <p>15 seconds</p>   <p>Window turns from clear to yellow when the injection is done.</p> <p>"click"</p> 